

# OYSTA

## OystaRio



- Fall Detection
- Safe Zones
- Auto Answer
- SOS Button
- Waterproof
- Audio Reminders
- Easy Charging
- Location Updates
- Two-way Voice Calls

### Rio Highlights

The Oysta Rio gives Vulnerable Independent People (VIPs) a discreet safety net so they feel confident to leave their home.

The Oysta Rio can be carried in a pocket, on a keyring, or on a lanyard. It has an SOS button, that, when pressed will notify a chosen contact in the care circle or Alarm Receiving Centre (ARC.)

The Oysta Rio connects to our care platform, IntelliCare™, so the care circle will be able to set safety parameters, locate VIPs and if needed, direct them home or direct help to their location.

It provides peace of mind to VIPs that they can remain independent and safe, and their care circle stays informed and available at a push of a button.

### The Oysta Rio is ideal for:

- Care Homes who would like to easily locate their care residents within the care home and if they are outside the care home.
- VIPs that like to be active but have started to experience memory loss/ disorientation and may need locating.
- People with learning difficulties who would benefit from the additional safety net whilst they gain confidence in their independence.
- VIPs who want the comfort of a telecare device but want it to be discreet.

### Benefits to Vulnerable Independent People (VIPs)

**One Touch SOS Button:** An obvious easy to use button for when in distress or simply needing peace of mind. The Oysta Rio connects the VIP to their chosen emergency contact or Alarm Receiving Centre (ARC) and makes a hands-free, two-way voice call.

**In the event of a fall,** an SOS call automatically be activated.

**Two-way Voice Calls:** The Rio has a loudspeaker and sensitive microphone to allow clear voice communication.

**Audio Reminders** for when the device has low battery to prompt the VIP to charge the device.

**Cradle Charging** means that the VIP only needs to drop it into the charging cradle to charge the device. The cradle also has an additional SOS button.

**Waterproof:** Should the VIP have trouble in the shower, they can press the SOS button for assistance.

**Auto Answer Option** automatically answers calls to the Oysta Rio. No need to press a button, saving time and effort in an emergency or a welfare check situation.



## Benefits to the Care Circle

### Locator Beacons

With Oysta Rio's Bluetooth Locator Beacons, you can identify via IntelliCare™ your VIP's location within a metre and whether they are at home or away from home.

This is particularly useful for care homes. Should residents leave the care home, you'll be notified when they leave the building and when they return.

### Fall Detection

Should the VIP have a fall, it will automatically connect them to their chosen contact or Alarm Receiving Centre (ARC)

### Safe Zones

Safe Zones are areas where the VIP is likely to travel. Should the VIP travel outside those safe zones, anytime or during specific times of the day or night, a care alert is triggered.

### WeatherWatch

A significant step towards preventative care; WeatherWatch notifies the care circle if the VIP is outside after dark or in extreme weather conditions.

You can conduct a welfare check before problems occur, and if needed, guide them home using the location feature on IntelliCare.

### Whitelist of Contacts

Only authorised callers from the care circle can contact the Oysta Rio which eliminates spam calls to the VIP. The whitelist can be changed at any time to add/remove contacts.

### Status Alerts

Through IntelliCare, you'll be notified if your VIP's device is switched on/off, has a low battery or if there has been 'no activity', indicating that the device has not been worn or used.

## OystaCare button

Should your VIP need assistance, you can guide help directly to their location. The OystaCare button can be worn discreetly. Should the VIP need assistance, they can press the button and through the Rio, it will alert the care circle to send assistance.



## Intelligent Care Technology

Our devices and sensors link into our easy-to-use care platform, IntelliCare. From IntelliCare you gain clear insight into your VIPs activity and patterns and receive notifications when activity falls outside of normal parameters. This will enable a clearer picture of your VIPs wellbeing and ensure you are best positioned to provide proactive and reactive care.

| SPECS              |   |
|--------------------|---|
| Dimensions         | 62mm*47.9mm*18mm  |
| Weight             | 60g   |
| FEATURES           |   |
| Microphone/speaker | Built-in Microphone & Speaker                               |
| Signal             | 4G LTE Supported Bands: Europe, Australia and North America |
| GSM Freq.          | Quad-band 850/900/1800/1900 MHz                             |
| GPS chipset        | Support GPS and GLONASS / AGPS                              |
| Backup Battery     | Rechargeable, 3.7V,1000mAh                                  |
| Battery life       | Up to 72 to 240 hours under normal usage.                   |
| Charging voltage:  | 5V DC   |
| Operating temp     | -20°C to +80°C for working -30°C to +70°C for storage       |
| Waterproof         | Follow the IP67 Standard                                    |

Call us on  
**01295 530 101**  
to book a live  
demonstration  
of the Oysta Rio  
and IntelliCare™



# OYSTA

Version: Sep2022

**OYSTA TECHNOLOGY**

Home Ground Barn, Pury Hill Business Park, Towcester, NN12 7LS, UK.

Tel +44 (0) 1295 530 101 Email [info@oysta-technology.com](mailto:info@oysta-technology.com) [www.oysta-technology.com](http://www.oysta-technology.com)

